

SPEED

RELIABILITY

SPECIFICITY



MICA
Advance

E. coli




Escherichia coli

Drinking water, including both bottles and tap water

- ✓ **Detects exclusively *Escherichia coli* (*E.coli*)**
- ✓ **Procedure based on the standard method, ISO 9308-01 with automatic detection**
- ✓ **6 hours of incubation only**, compared to 24 hours with the standard method
- ✓ **Up to 500 analyses per 8 hour period**
- ✓ **No additional confirmation required**
- ✓ **Limit of detection: 1 CFU* / filtered volume**

CE | UK
CA

Protocol based on the **ISO 9308-01** method



*CFU: Colony Forming Unit

More info

> diamidex.com